

Easy Dider I can go out by m

With the Easy Rider I can go out by myself and enjoy cycling again. The seat is very comfortable.

David





Easy Rider

Tricycle

The Easy Rider is a tricycle for adults with optimal riding characteristics that make for a pleasant cycling experience. Due to the unique frame with a low step-through, the bike has a low centre of gravity. This provides a safe and stable feeling.

This tricycle has a comfortable and easy to adjust seat. The forward pedalling motion and the suspension in the frame give ultimate comfort while cycling. The Easy Rider is therefore suitable for both long, and short trips.

The modern and aerodynamic design with a hydroformed frame gives the Easy Rider its sporty appearance. With the optional pedal assist (smart e-bike), cycling becomes even easier!

Benefits

- ✓ Suspended frame for optimal comfort
- ✓ Very easily manoeuvrable due to its small turning radius
- $\checkmark \;\;$ Comfortable seat with tailbone relief and adjustable backrest
- ✓ Suitable for up to 140 kg (reinforced frame up to 180 kg)
- ✓ The seat and handlebars are easily adjustable

Popular options/accessories

- A variety of frame sizes (M and L)
- Different suspension strengths for optimal comfort
- Stick holder or basket for carrying essentials
- A seat belt, mirror, armrests, or foot fixation for extra security
- Extra battery pack or battery with a higher capacity for an increased range

