## **FunTrain**

### Side-by-side tandem trailer

Cycling together becomes even more special with the FunTrain. With this side-by-side tandem trailer, which you attach to a Fun2Go using a towbar, one main driver can cycle with three passengers.

The main driver sits on the duo bike and has full control over the FunTrain since they are the only one who can steer and brake. The trailer has an automatic overrun brake, which ensures safe and efficient braking. On the FunTrain, you sit next to each other, making it easy to communicate and enjoy the scenery.

The retractable nose wheel allows for easy transportation of the detached FunTrain.







# It is nice to be able to go out with three residents. They always look forward to the rides we make.

#### Benefits

Kim

- ✓ One supervisor can cycle with three passengers
- ✓ Easy to connect to a Fun2Go duo bike
- ✓ Good visibility and easy communication with side-by-side seating
- Comfortable seats with backrest that are easy to adjust
- $\checkmark \;\;$  Disconnected FunTrain is easily movable with a retractable nose wheel

### Popular options/accessories

- Retractable footrest for passengers
- Seatbelts and armrests for extra safety
- Pedal adjustment, such as foot fixation
- Bike cover for covering the FunTrain
- Walking stick holder

