

EASY RIDER

Comfortable tricycle with seating for adults



The Easy Rider tricycle for adults is clearly different from other tricycles. It has a unique frame with low step through and a stable lower sitting position. The low centre of gravity, seating with a back support and steering in front, give the rider a safe and secure feeling. The seat height and seat angle are adjustable. The forward pedaling motion and the frame suspension ensure a perfectly comfortable ride.

The special cycling position lays hardly any burden on the back, neck, shoulders and wrists. Virtually anyone can ride this cycle without getting habituated to it. The Easy Rider is available with various options such as single hand brakes, a walker holder, foot fixation and an electric engine for power assist.

SUITABLE FOR:

Adults who want to ride comfortably and without fear on a beautifully designed tricycle.

SPECIFICATIES

Inside leg length	69 - 95 cm
Step through height	38 cm
Cycle weight	+/- 43 kg
Cycle length	210 cm
Cycle width	75 cm



Seat



Suspension



Divisible
frame

