

Technische Spezifikationen

	Easy Rider (M)	Easy Rider Compact	Easy Rider Compact Small	Easy Sport	Easy Sport Small	Maxi	Maxi Comfort	Midi	Mini	Husky	Viktoria	Viktor	Easy Go	Balance	City	Fun2Go	FunTrain	Twinny	Twinny Plus	Kivo	Kivo Plus	OPair	VeloPlus	Chat
Rahmenhöhe (cm)	-	-	-	-	-	53 ¹	51 ¹	41 ¹	30 ¹	23 ¹	41	50	-	51 ¹	47 ¹	-	-	52 ¹	52 ¹	30 ¹ _V 52 ¹ _H	30 ¹ _V 52 ¹ _H	52	52	52
Einstiegshöhe (cm)	35	32	32	52	41	28	23	25	20	12	18	16	29	22	22	-	-	37 _V 44 _H	37,5 _V 45 _H	36 _V 43 _H	36 _V 45 _H	45	28	45
Gewicht (± kg) <small>(Ohne Elektrosystem, außer Chat und Easy Go)</small>	49	39	37,5	39	36	29	32	27	21	11	30	30	64	26	9	76	76	35	44	39	48	58	64	97
Max. Gewicht Gepäckssystem (kg)	20	10	10	20	20	20	20	20	20	20	15	15	20	23	-	20	20	20	20	20	20	20	20	20
Länge (cm)	225 ¹ _{max.}	173 _{max.}	173 _{max.}	225 _{max.}	186	180	194	168	136	98	155	173	166	182	113	199	193	259	261	244	246	214	253	232
Breite (cm)	78	73	73	80	77	75	75	75	60	59	80	89	72	66	55	113	113	66	75	66	75	76	110	111
Option: Teilbarer Rahmen	✓	-	-	-	-	✓	-	✓	-	-	-	-	-	-	-	-	-	✓	✓	-	-	✓	✓	-
Max. Benutzergewicht (kg) (pro Person)	140 ²	100	80	120	80	120	120	90	65	65	110	110	140	120	100 ²	120	120	120	120	120	120	120	110/140 ³	120/200 ³
Innenbeinlänge (cm)	65-88 ¹	63-81	52-72	65-105 ¹	47-73 ¹	80-95 ¹	76-91 ¹	66-80 ¹	47-59 ¹	35-47 ¹	63-80	75-93	78-90 ¹	76-91 ¹	70-95 ¹	67-93 ¹	67-93 ¹	77-89 ¹	77-89 ¹	50-62 ¹ _V 77-89 ¹ _H	50-62 ¹ _V 77-89 ¹ _H	80-95	80-95	80-95

¹ = Anpassungsfähig ² = Verstärkte Ausführung möglich

³ = VeloPlus: Fahrer max. 110 kg, Passagier + Rollstuhl max. 140 kg* (*Verstärkter Frontbereich max. 200 kg).
Chat: Fahrer max. 120 kg und Passagiergewicht im Sitz max. 200 kg.

V = Vorne
H = Hinten