

FUNTRAIN

Double ride trailer



The Funtrain is an automatic, hydraulic braked double rider cycle trailer, which can be attached behind an existing or new Fun2Go side by side bike. In this manner, an attendant can safely ride with three co-riders or usually with one co-rider on the Fun2Go itself.

A tow bar is simply fitted on the existing Fun2Go. The Funtrain has a safety lock which cannot come off on its own and an automatic, hydraulic, overrun brake. The driver of the front double rider bike thus has 100% control over the Funtrain and Fun2Go.

The driver is the only one who steers and brakes, all the four riders can pedal along. Optionally, it can also be determined for the co-riders individually whether they must pedal, with the same speed or different speed or need not pedal at all. The Funtrain itself has a foldable support wheel to move the detachable double rider bike trailer easily.

You can even attach two Funtrains behind a Fun2Go so that six people can cycle together.

SUITABLE FOR:

Institutions, homes, and communities where one attendant wants to cycle with many individuals.

TECHNICAL SPECIFICATIONS

Seat height	61 cm
Step through height	none
Cycle weight	+/- 67 kg
Cycle length	197 cm
Cycle width	113 cm



Seatpillow
5 cm thick



Switchable
freewheel



Footfixation

