

# HUSKY

## Tricycle for young children



The Husky tricycle is specially designed for young children, who have difficulties with stability. This cheerful, light-running tricycle has a specially designed bottom bracket, making it possible for a child to pedal downwards. The crank normally has integrated crank-shorteners.

It is not possible to use an electric engine for power assist on the Husky.

### Standard specifications Husky

- Plastic chain guard
- Environment-friendly powder-coating
- Parking brake
- Steering limiter
- Five-year guarantee on frame
- Bell

### SUITABLE FOR:

Children aged 2 years and above.

### TECHNICAL SPECIFICATIONS

Inside leg length	35 - 47 cm
Frame height	23 cm
Step through height	12 cm
Cycle weight	+/- 11 kg
Cycle length	98 cm
Cycle width	59 cm



Push bar, back fixation and belt

