

# EASY GO

## Compact electric mobility scooter bike



The compact electric scooter bike is designed for people with stability problems or variable energy level. There are three options available on this tricycle:

- You can ride the trike by just pedaling.
- You can ride it by pedalling with electric pedal support.
- You can use 100% electric power assist.


When you fold the footrest down, you do not pedal anymore and use the cycle as a scooter. This can have variable speed levels. You can easily decide which type of support you want during the cycle ride. This cycle is quite compact and maneuverable and hence, also suitable to use for shopping. The adjustable double suspension gives you a stable and comfortable driving experience. There is no driving licence or any other driving ability certificate required to ride the Easy Go. However, there must be an insurance done for it.

### SUITABLE FOR:

Adults with a variable energy level, who want to cycle, but would also like to have the possibility of completely switching to electric scooter if they cannot cycle (any more).

### TECHNICAL SPECIFICATIONS

Inside leg length	75-90 cm
Step through height	28,5 cm
Cycle weight	+/- 64 kg
Cycle length	166 cm
Cycle width	71,5 cm



Adjustable suspension



Footrest



Throttle

