



# Maxi Comfort

## Tricycle with low step-through

Van Raam has an alternative version of the traditional tricycle. The Maxi Comfort tricycle has an extra low and wide step-through. This makes getting on and off the bike very easy.

On this tricycle, the seat is not placed directly above the bottom bracket, but slightly behind it. This allows for optimal movement of the knee while sitting slightly lower. Because the centre of gravity is also a bit lower, the tricycle is very comfortable.

Are you looking for a two-wheeled bike with a low step-through? The Balance might be a good option.



*"The Maxi Comfort tricycle gives me the security and confidence to go out." - Gerda*

### Benefits

- ✓ Easy to get on and off due to the low step-through
- ✓ Ergonomic and comfortable position
- ✓ More power due to forward pedalling motion
- ✓ Medium-high seating position for a good overview
- ✓ Lighter steering compared to a standard tricycle

### Popular options/accessories

- Pedal assist (Smart E-Bike)
- A variety of frame sizes
- Backrest (with belt) or foot fixation for extra security
- A basket or walking stick holder to bring along your necessities
- 8-gear hub with freewheel or backpedal braking

