

# CITY

## Walking support



The City walking support is suitable for people who have difficulty in walking or prolonged standing. A wheelchair is not yet required and you do not prefer a walker. Your body weight is borne by the seat of the walking support, thus relieving your leg joints. Walking takes less effort, hence you can travel over longer distances. Besides, the walking support is foldable and light-weight, so you can easily carry it along in a car, bus, train or plane.

### SUITABLE FOR

People having difficulty in walking or prolonged standing Walking support.

### SPECIFICATIES

Inside leg length	70 - 95 cm
Frame height	47 cm
Step through height	22 cm
Cycle weight	+/- 9 kg
Cycle length	113 cm
Cycle width	55 cm



Small frame

