

Technical specifications

	Easy Rider M	Easy Rider L	Easy Rider Junior	Easy Sport	Maxi	Maxi Comfort	Midi	Mini	Husky	Viktoria	Viktor	Easy Go	Balance	City	Fun2Go	FunTrain	Twinny	Twinny Plus	Kivo	Kivo Plus	OPair	VeloPlus	Chat
Frame height (cm)	-	-	-	-	53 ¹	51 ¹	41 ¹	30 ¹	23 ¹	42	50	-	51 ¹	47 ¹	-	-	52 ¹	52 ¹	30 ¹ front 52 ¹ rear	30 ¹ front 52 ¹ rear	52	52	52
Entry height (cm)	35	35	35	31	28	23	25	20	12	19	16,5	28,5	23	22	-	-	37 front 44 rear	37,5 front 45 rear	36 front 43 rear	36 front 45 rear	45	28	45
Weight of the bike (± kg) <small>(Pedal assist excluded)</small>	49	51	36	39	29	32	27	21	11	30	30	64	22	9	69	67	35	44	39	48	58	64	97
Max. weight luggage carrier (kg)	20	20	20	20	20	20	20	20	20	15	15	20	20	-	20	20	20	20	20	20	20	20	20
Length of the bike (cm)	225 max.	242 max.	175	195-225	180	194	168	136	98	155	173	166	182	113	200	197	259	261	244	246	214	253	232
Width of the bike (cm)	78	78	75	80	75	75	75	60	59	80	89	71,5	66	55	113	113	66	75	66	75	76	253	111
Divisible frame available	✓	✓	✓	-	✓	-	✓	-	-	-	-	-	-	-	-	-	✓	✓	-	-	✓	✓	-
Max. user weight (kg)	140 ²	140 ²	80	120	120	120	90	65	65	110	110	140	120	100 ²	120 ²	120	120	120	120	120	120 ³	110/140 ³	120/200 ³
Inner leg length (cm)	65-88	85-108	47-73 ¹	65-105 ¹	80-95 ¹	76-91 ¹	66-80 ¹	47-59 ¹	35-47 ¹	63-80	75-93	78-90 ¹	76-91 ¹	70-95 ¹	75-90 ¹	75-90 ¹	77-89 ¹	77-89 ¹	50-62 ¹ front 77-89 ¹ rear	50-62 ¹ front 77-89 ¹ rear	80-95	80-95	80-95

¹ = Customisable ² = Reinforced version available

³ = OPair: Both the driver and passenger max. 120 kg. VeloPlus: Driver max. 110 kg. Passenger + wheelchair max. 140 kg* (* Reinforced front max 200 kg). Chat: driver max. 120 kg and total passenger weight max. 200 kg.